Cuff Type	Corresponding disinfectant	Disinfection method	Result
SlimFit Cuff • CUF-F-LA • CUF-F-SL • CUF-F-A • CUF-F-SA Resin with D ring	Chlorhexidine gluconate aqueous solution (0.5%) Benzalkonium chloride aqueous solution (0.1%) Neutral detergent (500-fold dilution)	Hold the gauze soaked in the chemical in your hand, apply a load of 1 kg, and wipe the target area 550 times.	No effect to the cuff is confirmed.